

Tual[®] unguento

Tual[®] unguento è indicato in numerose condizioni di alterata reattività cutanea.

Per le sue caratteristiche antimicrobiche, disinfettanti, cicatrizzanti, antinfiammatorie ed emollienti Tual[®] è utile in caso di:

- Suture • Dermatiti • Irritazioni cutanee • Follicoliti • Cheiliti
- Post peeling • Secchezza della mucosa nasale (croste nel naso)

MODO D'USO:

Applicare la quantità necessaria di unguento ripetendo più volte al giorno, se necessario, fino a completa normalizzazione. Al tubetto può essere avvitata la cannula in dotazione per un'applicazione precisa ed efficace anche in punti difficili (es. cuoio capelluto, mucose dei genitali esterni). Evitare il contatto con gli occhi. Richiudere con cura il tubetto immediatamente dopo l'uso.

Uso esterno, cosmetico, conservare in luogo fresco e asciutto.

Tual[®] unguento, due to its components, is indicated in several conditions of impaired cutaneous reactivity.

USE:

Apply the necessary amount of ointment, if necessary repeat several times a day, until complete normalization of the skin. The tube can be equipped with the supplied cannula for a precise application even in difficult areas (for example scalp or mucosa of the external genitals). Avoid contact with eyes. Carefully close the tube immediately after use.

External use, cosmetic, store in a cool dry place.



Via S. D'Orsenigo, 5 - 20135 Milano
www.gamfarma.it